

CAUTION

- Instrumentation of the cervical os may result in vasovagal reactions, including fainting. Have the patient remain supine until she feels well, and have her get up with caution.
- If there is clinical concern, exceptional pain or bleeding during or after insertion, appropriate steps (such as physical examination and ultrasound) should be taken immediately to exclude perforation.
- The device is for single use only. Reusing a device might cause infections and risk of misplacement (including perforation) and should be avoided.

Instructions to be given to the patient by the doctor

- Invite the patient for a visit after her first menses to confirm that the IUB™ is still in place.
- If she feels any of the following symptoms she should contact her healthcare provider:
abdominal or pelvic pain, cramping or tenderness; unusual or malodorous discharge; unexplained vaginal bleeding; unexplained fever, chills; painful sex; a missed period; feeling that the length of the threads changed or feeling any other part of the IUB™ besides the threads; an allergic reaction.
- Patient should contact her healthcare provider if:
- Patient thinks she is pregnant; becomes HIV positive or her partner becomes HIV positive; might be exposed to sexually transmitted diseases (STDs).
- If partial expulsion or expulsion occurred, the patient is not protected from pregnancy.


Patient Follow-up

- Following placement, examine the patient after her first menses to confirm that the IUB™ is still in place. You should be able to see or feel only the threads. If the IUB™ has been partially or completely expelled, remove it. You can place a new IUB™ if the patient desires and if she is not pregnant. Do not reinsert a used IUB™.
- Evaluate the patient promptly if she complains of any of the following: abdominal or pelvic pain, cramping, or tenderness; malodorous discharge, bleeding, fever or a missed period.
- The length of the visible threads may change with time. However, no action is needed unless you suspect partial expulsion, perforation, or pregnancy.
- If you cannot find the threads in the vagina, check that the IUB™ is still in the uterus. The threads can retract into the uterus or break, or the IUB™ may have perforated the uterus or expelled. Radiography or sonography may be required to locate the IUB™. If there is evidence of partial expulsion, perforation or breakage, remove the IUB™.

How to Remove the IUB™

The IUB™ should not remain in the uterus for more than 5 years.

- Prepare sterile gloves and sterile forceps. Remove the IUB™ with forceps, pulling gently on the exposed threads.
- If the threads are not visible, determine location of the IUB™ by ultrasound.
- Inspect to assure the integrity of the IUB™, specifically to the presence of the leading and trailing copper balls (see image below).



- In case of absence of visible threads or breakage of the IUB™ removal can be difficult. Analgesia and cervical dilation may assist in removing the IUB™. An alligator forceps or other grasping instrument may be helpful. Hysteroscopy may also be helpful.

You may immediately insert a new IUB™ if the patient requests so and has no contraindications.

How is the IUB™ Supplied

The IUB™ is available in cartons of 1 (one) sterile unit. Each IUB™ is packaged in a sterile pouch together with an insertion tube, a slider and a push rod.

The IUB™ is supplied sterile. Method of sterilization is ethylene oxide.

Shipping and Storage Conditions

Store the IUB™ packaging in a dry environment at 15°C to 30°C. In these conditions the IUB™'s shelf life is 3 years.

Short term transportation of the IUB™ packaging should be limited to a temperature between -18°C and 55°C.

- SINGLE USE, DO NOT RESTERILIZE.
- NEVER RE-INSERT A USED IUB™.
- NEVER USE AN IUB™ IF THE PACKAGE IS DAMAGED OR OPEN.
- DO NOT USE PAST THE EXPIRY DATE.
- DISPOSE OF USED IUB™ AND ITS COMPONENTS USING BIO-HAZARD DISPOSAL PRACTICES.

	Catalogue number
	Batch code
	Date of manufacture
	Use-by date
	Manufacturer
	Authorized representative in the European Community
	Sterilized using ethylene oxide
	European Conformity mark Notified Body: LNE-G-MED (0459)
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged
	Temperature limit
	Keep dry
	Consult instructions for use
	Caution
	MR Conditional
	Recycle package after use

MRI Safety Information

Non-clinical testing demonstrated that the IUB™ is MR Conditional.

A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.
- Under the scan conditions defined, the IUB™ is expected to produce a maximum temperature rise of 1.4 °C after 15-minutes of continuous scanning.

Artifact Information

In non-clinical testing, the image artifact caused by the IUB™ extends approximately 2 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.



IUB™ is a trademark of OCON Medical Ltd.
Ballerine is a registered trademark of OCON Medical Ltd.
www.oconmed.com

European Representative:
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Manufacturer:
OCON Medical Ltd.
Hamikso'ot Blvd. 14
Modi'in 7178095, Israel
Tel: +972 72 21 50 105

Imported and distributed by:

Consilient Health Limited
5th Floor, Beaux Lane House,
Mercer Street Lower, Dublin 2,
Ireland

Initial date of authorization: December 2014